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Deparaffinization

The following steps describe the deparaffinization procedure for a single maximum 10 µm thick section of formalin-fixed, paraffin-embedded tissue in a 1.5 ml reaction tube.

- 1.To one maximum 10 µm section in a 1.5 ml reaction tube, add 800 xylene, and vortex briefly in several intervals.
- 2.Add 400 µl absolute ethanol and vortex briefly. Centrifuge for 2 minutes at maximum speed (16,000 × g) and discard supernatant. Be careful not to dislodge the pellet.
- 3.Add 1 ml absolute ethanol and vortex briefly. Centrifuge for 2 minutes at maximum speed (16,000 × g) and discard supernatant. Be careful not to dislodge the pellet.
- 4.Blot the tube briefly onto a paper towel to get rid of ethanol residues. Dry the tissue pellet for 10 minutes at +55°C (open tubes). If required, dry for up to a maximum of 20 minutes. Proceed to Step 1 of the RNA isolation protocol

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Purity

Purified RNA is free of DNA, nucleases, and all cellular and sample components that interfere with RT-PCR, according to the current Quality Control Procedures.

The purity is determined by calculation of OD 260/280nm .Ratio 260/280 Eluate: 1.8 - 2.0 RNA Integrity and size distribution

The size of the isolated RNA fragments is analyzed using the Agilent Bioanalyzer. Ninety percent of RNA fragments are greater than 100 nucleotides (exemplary data). The RIN number is usually very low for FFPE-derived RNA samples due to fragmentation. RIN should be >1.4 as a general rule, but even RNA samples with no RIN may work if your PCR amplicons are small.

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Test Principle

FFPE tissue samples are disrupted and homogenized during incubation with RNA Tissue Lysis Buffer and Proteinase K. Nucleic Acids (NA) bind in the presence of a chaotropic salt specifically to the surface of glass fibers pre-packed in the High Pure Purification Filter Tube.

RNA is purified in a series of rapid "wash-and-spin" steps to remove salts, proteins, and cellular component Residual DNA is digested on the column by adding DNase I.

Finally, a low-salt elution releases the RNA from the glass fiber. The process does not require RNA precipitation or organic solvent extraction, ideal for rapidly purifying many samples simultaneously

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